



IMCO Instrument Lubricant Bath



Non-silicone and non-toxic. Antimicrobial, water soluble, steam permeable. Rust inhibiting.

1080-T2-IMC Instrument Bath, 4 gal btls/cs

IMCO Enzy-Clean Enzymatic Detergent



Enzy-Clean Enzymatic Detergent disintegrates blood products.

Neutral pH, bacteriostatic, double enzyme. Non-corrosive, safe for scopes. Cleans hard to reach surfaces. Can be used in ultrasonic cleaners.

1080-B9-IMC Enzy-Clean, 4 gal btls/cs

IMCO Instrument Detergent



Biodegradable with neutral pH. Hemolytic to remove blood. Rust inhibitor, rinsing agent. Highly concentrated. Four gallon bottles per case.

1080-B2-IMC High Suds Detergent for Manual Cleaning

1080-B1-IMC Low Suds Detergent for Ultrasonic Cleaners

1080-B3-IMC Powdered Cleaner, 5 lb., 6/cs



Medica 28 Plus

Medica 28 Plus is a new formulation with accelerated 25-minute kill time for faster high-level disinfection of instruments. Medica 28 Plus uses a lower glutaraldehyde level (3.0%) than other extra strength products to reduce odor and vapor level in the operator. Medica 28 Plus has a minimal surfactant level to protect lensed instruments including endoscopes and a new solid activator to assure a mild pleasant smell and to maintain a consistent efficacy level through 28 days of continuous use. Contains rust and odor inhibitors.

2459-A007-02 4 gallons per case

IMECO
Private Label
Low cost without compromising quality

**Saving money
is no small task -
Let us help!**

Look for the IMCO logo
throughout this catalog.

www.imcoinc.com

The Products Listed In This Catalogue Represent Only A Portion Of Those Available. Please Call Our Office Or See Your Sales Representative.

MetriZyme® Dual Enzymatic Detergent

- Safe for instruments
- Three times more proteolytic enzyme activity than the competition
- Less corrosive than distilled water
- Safe for users

4397-10-4000	MetriZyme® gallon, 2 bottles/case
4397-10-4005	MetriZyme® quart, 4 bottles/case
4397-10-4010	MetriZyme® 1/2 gallon, 4 bottles/case
4397-10-4015	MetriZyme® 15 gallon, 1 drum/case
4397-10-4030	MetriZyme® 30 gallon, 1 drum/case
4397-10-4050	MetriZyme® 5 gallon, 1 drum/case
4397-10-4055	MetriZyme® 55 gallon, 1 drum/case
4397-46-1000	Dispensing Pump – 1 gallon, each
4397-46-1010	Dispensing Pump – 1/2 gallon, each
4397-46-1020	Dispensing Pump – quart, each

MetriSponge® Sponge Saturated with MetriZyme®

- Contoured sponge simplifies the effective pre-cleaning of surgical instruments and scopes
- Unique shape enables thorough cleaning of cylindrical surfaces without scratching

4397-10-4025	MetriSponge® 25 each/box, 4 boxes/case
--------------	--

DetergeZyme® Enzymatic Presoak and Cleaning Solution

- Non corrosive
- Easy-to-use
- Safe for automated systems

4397-46-1000	Dispensing Pump – 1 gallon, each
--------------	----------------------------------

EmPower Foam Foaming Enzymatic Spray

- EmPower Foam covers the instruments, reducing splashing and messy spills
- Ready-to-use in a convenient sprayer
- Breaks down bio-burden on instruments
- Speeds up the cleaning process

4397-10-4224	EmPower Foam, 24 ounce spray, 12 bottles/case
--------------	---



EmPower Dual Enzymatic Detergent

- EmPower works fast to remove blood, protein and other proteinaceous debris
- Economical concentrated product
- Safe for instruments
- Less corrosive than distilled water
- Contains two gold standard proteolytic enzymes for more effective cleaning
- Non-foaming
- Dilution 1:128



4397-10-4100	EmPower, gallon, 4 bottles/case
4397-10-4115	EmPower, 15 gallon, 1 drum
4397-10-4130	EmPower, 30 gallon, 1 drum
4397-10-4150	EmPower, 5 gallon, 1 bottle/case
4397-46-1000	Dispensing Pump, 1 gallon, each

MetriWash™ Instrument Detergent Concentrate

- Cost effective
- Improves cleaning and rinses freely
- Non-corrosive
- No phosphates, chlorides, or hydroxides
- Neutral pH

4397-10-3300	MetriWash™ gallon, 4 bottles/case
4397-10-3350	MetriWash™ 5 gallon 1 drum/case
4397-10-3315	MetriWash™ 15 gallon, 1 drum/case
4397-10-3330	MetriWash™ 30 gallon 1 drum/case
4397-46-1000	Dispensing Pump – 1 gallon, each

MetriClean®2 Low Foam Instrument Cleaner & Polisher

- Low foaming
- Cleans blood and inorganic soils
- Removes Stains
- Increases instrument functionality
- Alkaline pH
- For use in all kinds of water conditions
- Restores shine to surgical instruments

4397-10-8100	MetriClean®2 gallon, 4 bottles/case
4397-10-8115	MetriClean®2 15 gallon, 1 drum/case
4397-10-8130	MetriClean®2 30 gallon, 1 drum/case
4397-10-8150	MetriClean®2 5 gallon, 1 drum/case
4397-46-1000	Dispensing Pump – 1 gallon, each

Miltex Instrument Cleaners

- 4484-3-725 Surgical Instrument Cleaner, 1 gallon pump (3.78 liter) bottles (case of 4). Contents of 1-gallon container makes 512 gallons of solution. Neutral pH and phosphate-free liquid concentrate formula for manual and ultrasonic cleaning or soaking of surgical/dental instruments
- 4484-3-726 Surgical Instrument Cleaner, 1/4 oz. (7.4 ml) packets (dispenser box of 80). Contents of 1-box makes 80 gallons of solution. Neutral pH and phosphate-free liquid concentrate formula for manual and ultrasonic cleaning or soaking of surgical/dental instruments
- 4484-3-755 EZ-Zyme All Purpose Enzyme Cleaner, 1 gallon pump (3.78 liter) bottles (case of 4). Contents of 1-gallon container makes 170 gallons of solution. Neutral pH multi-enzyme liquid concentrate ideal for instrument presoak, ultrasonic cleaning, laundry presoak, and manual cleaning.
- 4484-3-760 Instrument Prep Enzyme Foam, 24 oz. (0.71 liter) spray bottles (case of 12). Ready to use foaming spray for pre-cleaning of soiled instruments and scopes (rigid & flexible). Ideal for applications where instrument cleaning is not immediately available.



MetriCide® High Level Disinfectant/Sterilant

- 2.6% Glutaraldehyde
 - Compatible with medical equipment and instrumentation
 - Contains no surfactants
 - No dilution required
 - Ideal for scopes
 - 14 day product
- 4397-10-1400 MetriCide® gallon, 4 bottles/case
 4397-10-1405 MetriCide® quart, 16 bottles/case
 4397-10-303 MetriTest™ Strips 1.5%, 60 strips, 2 bottles/case

MetriCide® 28 High Level Disinfectant/Sterilant



- Recommended for immersible medical and heat sensitive equipment and instruments
 - Contains a surfactant
 - Contains a rust-inhibitor to reduce corrosion
 - No dilution required
 - 28 day product
- 4397-10-2800 MetriCide® 28 gallon, 4 bottles/case
 4397-10-2805 MetriCide® 28 quart, 16 bottles/case
 4397-10-304 MetriTest™ Strips 1.8%, 60 strips, 2 bottles/case

MetriLube® Instrument Lubricant Concentrate

- Oil based lubricant
 - Extends life of instruments
 - Anti-corrosive
 - Non-greasy
 - Economical
- 4397-10-3400 MetriLube® gallon, 4 bottles/case
 4397-10-3450 MetriLube® 5 gallon, 1 drum/case
 4397-10-3415 MetriLube® 15 gallon, 1 drum/case
 4397-10-3425 MetriLube® 24 ounce spray, 12 bottles/case
 4397-10-3430 MetriLube® 30 gallon, 1 drum/case
 4397-46-1000 Dispensing Pump – 1 gallon, each

MetriCide Plus 30® High Level Disinfectant/Sterilant

- 3.4% Glutaraldehyde
 - Extra-strength solution
 - Recommended for heat sensitive equipment and instruments
 - No dilution required
 - Contains a surfactant
 - 28-day product
 - * Using MetriTest 1.8% Test Strips is recommended before each use.
- 4397-10-3200 MetriCide Plus, 30 gallon, 4 bottles/case
 4397-10-3205 MetriCide Plus, 30 quart, 16 bottles/case
 4397-10-304 MetriTest Strips, 1.8%, 60 strips, 2 bottles/case

